

Body Armor Pro Term Foot Stump Orthosis

DARCO's Body Armor Pro Term is a provisional orthosis for use after Chopart and Lisfranc amputations. Transmission of the load to the proximal foot is ensured by the wedge-shaped outer sole and the inner sole. The robust cross ankle shaft stabilizes both the talocrural and talocalcaneal joints.

Built-in adjustable air cushions allow this orthopedic device to adjust to fully fit and support the stump of the foot. The ventral support tongue provides compression in an extended area thanks to the variable closure system. The mounted orthosis can be customized to meet the individual needs of the patient.

The two-piece insole, composed of an EVA rope and an EVA / PE sole, is interchangeable or customizable as required.

Product Features & Benefits:

- Wedge-shaped outsole and insole in dorsal extension ensure rear weight-bearing.
- The stable, ankle-spanning shaft stabilises the upper & lower ankle and prevents deformation of the stump.
- The blue inlay with integrated pump compresses & this maintains the shape of the stump.
- Perfect fit through air compression. Easy inflation and deflation using an integrated system.
- This conservative treatment can also be manually adapted for customised support to meet individual patient needs.
- Designed to be worn for approximately 10 weeks following an amputation to allow the foot stump to completely heal. Thereafter, a conventional, custom-made shoe is worn.
- Universal fit Left/ Right in two sizes. Suitable for the right and left foot.
- Colour: Black
- Instructions for use included with each product, also available as a download from our website.

Indications for Use:

- Diabetic foot syndrome
- PAOD (peripheral arterial occlusive disease)
- Trauma
- Used post-operatively for treatment following Chopart or Lisfranc amputation

Cautions:

- The use of the brace should be supervised by a physician or health care specialist. In case of any abnormal sensation such as tingling, burning or numbness, pain, swelling or skin irritation, discontinue use and seek medical assistance.

Warranty:

- 12-month cover for defects in materials and workmanship in products purchased from Algeos only.

Care and Maintenance:

- Remove the liner and wash it by hand in cold water with a mild detergent. Do not use a washing machine or dryer. Place the liner flat to dry. Wash the plastic parts with a damp cloth and mild soap.
- Reusable if not used on open wounds.

End of Use of the product:

Metal and Plastic materials may be recycled. Follow local environmental guidelines for waste disposal.

Order Codes:

Optional: Replacement Liner (NV7584, NV7585)

Code	Size	Length
NV7582	Medium	16 cm
NV7583	Large	19 cm

